


THE ULTIMATE CHOICE
FOR PERIPHERAL ARCHITECTURE



zephyr[®]cs

Cobalt Chromium Covered Peripheral Stent

CoCr STENT DESIGNED FOR LARGE VESSEL SUPPORT

- Cobalt Chromium Stent Platform with ePTFE covering
- Zigzag hybrid stent design for easy navigation
- C & S Shape flexible poly links provides high radial strength and minimizes stent foreshortening
- Lower strut thickness provides minimal crimping profile
- Better Conformity & Stability
- Reduces the chances of Neointimal growth
- Long Term Patency



SPECIFICATIONS:

Stent Covering Material ePTFE Tubing	Stent Material Cobalt Chromium Alloy	Strut Thickness 0.009"	Stent Links C & S Polylinks
Expanded Stent Diameter Ø19 mm – Ø25 mm	Unexpanded Stent Diameter 5.0 mm	Stent Foreshortening <10 %	
Sheath Compatibility: 14F ± 1F <small>(subject to PTA used)</small>	Minimum Crimping Profile 12F ± 1F	Recommended Guidewire 0.035"	

SIZE MATRIX:

Configuration (Number of Crowns)	Diameter (mm)	Length (mm)								
		29	34	39	44	49	59	69	79	
12 Crowns	19	✓	✓	✓	✓	✓	✓	✓	✓	✓
	20	✓	✓	✓	✓	✓	✓	✓	✓	✓
	21	✓	✓	✓	✓	✓	✓	✓	✓	✓
	22	✓	✓	✓	✓	✓	✓	✓	✓	✓
	23	✓	✓	✓	✓	✓	✓	✓	✓	✓
	24	✓	✓	✓	✓	✓	✓	✓	✓	✓
	25	✓	✓	✓	✓	✓	✓	✓	✓	✓

ORDERING INFORMATION:

Product Codes	Length (mm)							
	29	34	39	44	49	59	69	79
	ZPCS29	ZPCS34	ZPCS39	ZPCS44	ZPCS49	ZPCS59	ZPCS69	ZPCS79

* SLTL recommends use of BIB Balloon, Z-Med Balloon, Maxi-LD Balloon or Cristal Balloon with Zephyr stent family

Corporate Office - Sahajanand Laser Technology Ltd.
E-30, Electronic Estate, G.I.D.C., Sector-26, Gandhinagar - 382 028, Gujarat, India.
Tel : +91-79-6170-1616, Fax : +91-79-2328-7470
Cell : +91-90-9904-0080
Email : mkt@SLTL.com
Log on to :
www.SLTLmedical.com | www.SLTL.com

Manufacturing Site - Sahajanand Laser Technology Ltd.
A-8, Electronic Estate, G.I.D.C., Sector-25, Gandhinagar - 382 016, Gujarat, India.
Cell : +91-90-9904-0090, Fax : +91-79-2328-7470